


## Disinfectant Cleaner

A concentrated, quaternary ammonium disinfectant cleaner that effectively cleans as it disinfects. Also inhibits growth of mold and mildew and controls related odors. Each 10-ml. cartridge yields one quart of spray disinfectant cleaner. EPA registered as a Broad Spectrum disinfectant cleaning product.

## Directions

- 1 Hold container with trigger sprayer attached 6–8 inches from surface and lightly spray surface. Wipe clean with towel or cloth.
- 2 Save bottle and sprayer for re-use.  
(For use as an AI cleaner in food processing applications, use of this product must be followed by a potable water rinse.)

**NOTE: Pre-test on painted surfaces in small, inconspicuous area for colorfastness.**



- 3) ATTACH SPRAYER TO BOTTLE & SHAKE GENTLY.
- 2) INSERT CARTRIDGE INTO BOTTLE NECK.
- 1) COLD WATER FILL LEVEL.

NOTE: DO NOT OVERFILL BOTTLE, CARTRIDGE WILL NOT EMPTY.



Ideal for athletic facilities and hospitals.

## Specs

Dilution: 1 cartridge makes 32oz  
pH Level: 6.2 - 7.5  
Color: Green  
Fragrance: Citrus

## Features

- Bloodborne Pathogen Compliant.
- Inhibits growth of mold & mildew & mitigates related odors.
- Broad Spectrum Disinfectant.

## Recommended For:

- Floors
- Walls
- Countertops
- Tile
- Doorknobs
- Fixtures
- Plastic Surfaces